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Instruction for use Kera[®]waxDisc Press

Milling wax for the CAD/CAM manufacturing of wax forms

Indication	Wax-ups of veneers, inlay, onlay, crowns and bridges, etc. Wax-frames for casting or pressing technology, e.g. with Eisenbacher dental alloys, press ceramic or thermoplastic.
Contra-indication	Not approved for the use in the human mouth.
Application	Perform a CAD design with appropriate software. Kera@waxDisc Press can be processed with all open CAD / CAM milling systems. Milled objects can be cut out of the blank using a scalpel or suitable tool. After final inspection on the master model, sprues can be designed and embedded with an investment material (e.g. Kera@Vest).
Technical data	Drop point: Approx. 80 ° C - 90 ° C Flash point: > 178 ° C Density: 1.02 g / cm ³ Color: vanilla Ash residue: Burns without residue
Form and size	ø98,3 mm x 20 mm with step
Contraindications and side effects	Allergies to the product or its components may occur in isolated cases. If an allergy is suspected, request a list of composition from Eisenbacher Dentalwaren ED GmbH. The product contains waxy acid mixture: esters of montanic acids, waxy acid mixture approx. C24-C34, paraffin, microwaves as well as confinement of polymer waxes.
Caution	If material comes into contact with the eyes, rinse it out with water for several minutes. If symptoms persist, consult an eye specialist.
Storage conditions	Protect from moisture and light.
Disposal instructions	Consult the material safety data sheets or national regulations for disposal.
Packaging unit	1 Disc 98,5 X 20 mm

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.