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# Instruction for use Kera®S-Powder

## Intended use

**Kera®S-Powder** is a beryllium- and nickel free cobalt based bonding alloy powder for the production of crowns and bridges by laser melting (SLM/ALM).  
**Only for professional user!**

## General guidelines for handling

This instruction for use includes important processing steps and recommendations for **Kera®S-Powder**. The alloy is solely intended for processing in laser sintering systems. Please refer to the operating instructions of the system manufacturer.

## Indication

Crowns and bridges, implant supported superstructures and bars, tertiary structures, model cast dentures.

## Dosage form

**Kera®S-Powder** is available with a corn size of 10 - 45 µm ± 5%.

## CAD

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm<sup>2</sup>). Sharp edges and undercuts should be avoided.

## Thermal treatment

After the laser sintering process, the building boards have to pass heat treatment to minimize stresses due to the laser sintering process. For this purpose, a suitable furnace with inert gas (argon) or vacuum function should be used. Stress relieving without inert gas atmosphere can optionally be performed. Please consider that a treatment without inert gas atmosphere can lead to an increased oxide formation.

Description	Temperature [°C]	Time [min]
Heat up 1 with inert gas within 60 min.	460	42
Heat up 2 with inert gas within 45 min.	760	~ 55
Cooling 1 with inert gas	Cooling to 600, then open door	
Cooling 2, w/o inert gas	300 – room temperature	

## Remove Parts from the build plate

After heat treatment and cooling the plate the restorations can be removed by band saw, rotary instruments or piers.

## Soldering / Laser welding

**Kera®S-Powder** can be soldered with all suitable solder. **Kera®S-Powder** parts should not be soldered with gold or palladium solder. **Kera®S-Powder** is also ideally suitable for laser-welding.

## Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C with vacuum (cleaning firing). The frame needs to be sandblasted with aluminium oxide 110µm and 3-4 bar to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

## Handling conditions / Safety

**Metal powder or dust may cause irritation by inhalation and in contact with skin. During the handling with the Kera®S-Powder as well as while grinding and sandblasting of the produced units it's recommended to consider an adequate extraction system, goggles, gloves, protective clothing and a respirator with fine particle filter (type FFP3 – DIN EN 149). After working with the powder or restorations clean hands thoroughly.**

## Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

## Single-use

Powder residuals should be not processed for further use as a medical device.

## Disposal Instructions

Consult the material safety data sheets or national regulations for disposal. Dispose of **Kera®S-Powder** residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling

## Storage conditions

Store powder in tightly closed containers and keep dry.

## Quantity

Please consider the label on the package.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

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