PRODUKTNAME

Kera® Soft-Disc

BEZEICHNUNG DESCRIPTION / DESCRIPTIF

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Instruction for use Kera®Soft-Disc

Intended use / Description

Kera®Soft-Disc is a biocompatible CoCr sintering blank for CAD/CAM processing. The milled framework must be sintered under argon gas. **The product may only be processed by trained qualified staff!**

Indication

Crown and bridge frameworks in the anterior and posterior region with max. 2 pontics, free-end bridges with max. one pontic, cone/telescopic crowns.

Processing

Preparation recommendation	Connector cross section, anterior	Connector cross section, posterior	Wall thickness	Margin thickness
step or pronounced chamfer	6 mm ²	9 mm ²	0,4 – 0,5 mm	0,2 mm

- For restorations > 3 units, a sinter support analogous to zirconium oxide must be considered.
- The shrinkage factor is indicated on the blank and must be set in the CAM software parameters.
- Zirconium oxide or carbide cutters with the milling strategy for zirconium oxide or sinter metal are suitable for milling.
- The milling process is carried out dry under aspiration or with water cooling.
- With water cooling, dry the milled parts for 30 min at 120°C before sintering.

Information for sintering process

- Before sintering, the objects must be cleaned of milling dust with a brush.
- Make sure that the objects are completely embedded in the sinterball bed (approx. 1 - 2 mm distance from the margin).
- Bridge pontics, bar parts, etc. must not be covered with sinter beads.
- The objects must be placed in such a way that sinter shrinkage is not impeded.
- If several restorations are sintered together, they must not touch each other.
- After cooling to < 50°C, the finished parts can be removed from the furnace.

Sintering parameter

Programmable sintering furnace with argon		Amann GIRRbach	
Heating rate	10°C / min	Heating rate	Standard
Final temperature	1280 °C	Final temperature	Standard
Holding time at final temperature	1 h	Holding time at final temperature	Standard
Cooling from end-temperature	Unregulated up to 300°C in argon atmosphere, from 300°C in air up to 50°C	Cooling from end-temperature	Standard

Soldering

Kera®Soft-Disc can be soldered with all suitable solder. **Kera®Soft-Disc** parts should not be soldered with gold or palladium solder.

Preparation before ceramic veneering

The frameworks are ground with cross-cut carbide cutters; ensure smooth transitions; avoid material overlaps. The minimum thickness of the caps should not be less than 0.3 mm. It is recommended that the frameworks are blasted with at least 110 µm aluminium oxide at 3-4 bar and cleaned (steam cleaner). An oxide firing can optionally be carried out at 980°C under vacuum with a holding time of one minute.

Kera®Soft-Disc can be veneered with standard dental ceramics for CoCr. Please observe the process steps of the respective ceramic manufacturers.

Handling conditions / Safety

Metal dust is harmful to health. Always wear your personal protective equipment when handling the product. For milling, cleaning, finishing and sandblasting use suction (HEPA H filter) and respirator with FFP3 filter - DIN EN 149, gloves and goggles. Do not use compressed air for cleaning.

Contraindications and side effects

Incompatibilities with CoCr alloys are extremely rare if these instructions for use are observed. In case of a proven allergy to a component of this alloy, it must not be used for safety reasons. In exceptional cases, electrochemically induced local irritations are described. When using different alloy groups in the patient's mouth, galvanic effects can occur. Please inform your dentist about contraindications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

Disposal Instructions

Dispose of **Kera®Soft-Disc** residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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