**PRODUKTNAME**  
PRODUCT NAME / NOM DU PRODUIT  
**BEZEICHNUNG**  
DESCRIPTION / DESCRIPTIF

### Kera®starPEEK

Biokompatibles, teilkristallines Hochleistungspolymer  
Biocompatible, semi-crystalline high performance polymer  
Polymère haute performance biocompatible, semi-cristallin

**FARBEN**  
COLOURS / COULEURS

TC (zahnfarben / tooth-colour / couleurs de dent)  
SW (weiß / star white / blanc)

**ABMESSUNG**  
DIMENSION / DIMENSIONS

Ø 98,3 mm mit Stufe / with step / avec épaulement:  
12 mm, 16 mm, 20 mm, 25 mm  
Ø 99,5 mm ohne Stufe / without step / sans épaulement:  
12 mm, 16 mm, 20 mm, 25 mm

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# Instruction for use Kera®starPEEK

## Intended use

Kera®starPEEK is a milling disk with filled Polyetheretherketone (PEEK) for the production of temporary and permanent dental restorations. **Only for professional user!**

**Not intended for the production of implants!**

## General guidelines for handling

This instruction for use includes important processing steps and recommendations for Kera®starPEEK.

## Indication

Full anatomical crowns and bridges (max. 2 pontics with sufficient dimensioned connectors (min. 10 - 12 mm<sup>2</sup>)  
Anatomical reduced crowns and bridges for the composite veneering with sufficient dimensioned connectors (min. 10 - 12 mm<sup>2</sup>)  
telescopic crowns and bridges, secondary structures, gingivaformer, abutments, splints, denture base, bar supported structures.

## Frame design

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with light curing composite. The wall thickness should not be less than 0.5 mm. Choose a sufficient connector dimension for the anterior and posterior region (10 - 12 mm<sup>2</sup>).

## Milling

Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine. Kera®starPEEK is a thermoplast material and should be processed with water cooling. For milling without water cooling a dust extraction system have to be installed. Use only milling tools which are approved for the processing of polymers and thermoplastics.

## Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

## Preparation of the surface before veneering with composite / resin

It's recommended prior the veneering to sandblast the surface with aluminium oxide (min. 50 µm, 3-4 bar) and clean by steam cleaner. It is mandatory to use an additional chemical conditioning e.g. visio.link, bredent, Senden for the frame.

## Fixation

For insertion of permanent restorations it's approved to use composite cements (z.B. Panavia SA Cement, Kuraray Medical Inc.). Consider a sufficient mechanical retention by sandblasting (min. 50 µm, 3-4 bar) and careful cleaning by steam cleaner. Use a chemical conditioning (e.g. visio.link, bredent, Senden) on the surface.

## Handling conditions / Safety

**Dust may cause irritation by inhalation and in contact with skin. During grinding and sandblasting it's recommended to consider an adequate extraction system and goggles.**

## Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with PEEK are extremely rare. In case of a proven allergy against an ingredient it should not be used for safety reasons. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

## Single-use

Used discs should be not processed for further use as a medical device.

## Disposal Instructions

Dispose of Kera®starPEEK residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

## Storage conditions

Store dry and protected against light.

## Quantity

Please consider the label on the package.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

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