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Instruction for use KERA® LabPutty 85

Temperature resistant kneadable silicone for dental laboratories, addition curing.

Indication For the quick production of models for denture repairs and creation of silicon keys for finalizing denture works.

Contra-indication Do not use **KERA®LabPutty 85** on patients.

Properties

Component A	White
Component B	Blue
Mixing ratio	1 : 1
Mixture time	30 s
Working time	1:30 min
Setting time	3:30 min
Hardness Shore A	85
Detail reproduction	20 µm
Shrinkage	< 0,1 %

Application

- Dispense both components AB in equal parts as required (1: 1).
- Mix both components manually about 30 s, until a homogeneous color is obtained.
- **KERA®LabPutty 85** should be processed in 1:30 min.

Handling conditions / Safety

- Close containers carefully after use.
- Keep caps / lids / spoon always separated.
- When processing the material wearing of vinyl or nitrile disposable gloves is recommended.
- Processing at room temperature 23 ° C.
- Do not use **KERA®LabPutty 85** if the expiry date (see label) has elapsed.

Contraindications and side effects

Allergies to the product or its components may occur in isolated cases. If an allergy is suspected, request a list of composition from Eisenbacher Dentalwaren ED GmbH.

Caution

If material comes into contact with the eyes, rinse it out with water for several minutes. If symptoms persist, consult an eye specialist.

Storage conditions

Do not store above 25 °C. Keep out of children's reach.

Disposal instructions

Consult the material safety data sheets or national regulations for disposal.

Quantity

KERA®LabPutty 85 A - 5 Kg / 1,5 Kg
KERA®LabPutty 85 B - 5 Kg / 1,5 Kg

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.