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GHS07

# Instruction for use Light-Curing Tray-Material

Light-curing base material for the fabrication of individual impression trays and bite rims based on acrylates.

- Indication**
- Individual impression trays
  - Base for denture try ins
  - Base for bite rims
  - Base for bite registrations

**Contra-indication** All application that are not listed in the point indication.

- Preparations**
- Sketch the form of the base plate on the model.
  - Block out undercuts on the model with wax. For an individual tray, block out the whole model (e.g. with plate wax).
  - Isolate the plaster model with an alginate based isolation.

- Application**
- Adapt the base plate on the blocked out plaster model and cut e.g. with a scalpel in the desired previously sketched form. **Tip: When making an individual tray, the remaining material can be used for the handle.**
  - Cure the adapted plate in a light curing unit.

**Polymerization times:**

**Process the curing in light curing units with a wave spectrum of 350-470 nm for 3- 5 min.**

- Gently peel the individual tray from the plaster model and polymerize when needed from the bottom again.
- Iso-propanol can be used to remove the inhibition layer.
- Finish with carbide cutters and sandpaper for smoothing the edges.

**Technical data** Flexural strength: > 130 MPa

Processing under daylight: ca. 20 min

- Handling conditions / Safety**
- When processing the material it's recommended wearing vinyl or nitrile disposable gloves.
  - Avoid skin and eye contact with the unpolymerized material.
  - When grinding please use a dust extractor.
  - Do not use Light-curing tray material after the expiry date (see label).

**Contraindications and side effects** Allergies to the product or its components may occur in isolated cases. If an allergy is suspected, request a list of composition from Eisenbacher Dentalwaren ED GmbH. The material contains acrylates and photoinitiators.

**Caution** If material comes into contact with the eyes, rinse it out with water for several minutes. If symptoms persist, consult an eye specialist.

**Storage conditions** Keep out of children's reach. Protect against light.

**Disposal instructions** Consult the material safety data sheets or national regulations for disposal.

**Packaging unit** Light protected box with 50 plates. Color: Pink

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

GHS07