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Instruction for use Kera®NH

Intended use

Kera®NH is a beryllium free nickel based bonding alloy for the production of crowns and bridges.

Kera®NH should not be used in patients with a known sensitization to nickel. **This product contains Nickel! Only for professional user!**

General guidelines for handling

This instruction for use includes important processing steps and recommendations for **Kera®NH**.

Indication

For the production of crowns and bridges for the ceramic veneering.

Modelation

The modelation should be done with wax that fire without leaving residues under consideration of the standard rules of designs for dental technicians. The frame has to be designed in an anatomical reduced form. The wall thickness should be at a minimum of 0,4 mm to secure the flow out of the melt. Consider a sufficient connector (6-9 mm²). Avoid sharp edges and undercuts.

Sprue design and investing

We recommend the design of the sprue with a bar. The horizontal sprue should have Ø 4-5mm, the sprue to the restoration should have Ø 3 mm. Single crowns should be directly connected with a sprue of Ø 4 mm with a length of 15-20mm. Connect the sprue on massive areas e.g. palatal and avoid the center of the muffle. For the investing we recommend the use of **Kera®-Vest** (a phosphate-bonded investment from Eisenbacher Dentalwaren ED GmbH). The pre-heating process of the muffle to an end temperature of 900°C has proved to be successful. For more information please read the instruction for use of **Kera®-Vest**.

Melting and casting

Kera®NH should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material.

Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and the cast shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1430°C.**

Devesting and cleaning

Let the muffle cool down to room temperature (approx. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 110µm of aluminium oxide with 3-4 bar, then clean with a steam cleaner.

Soldering / Laser welding

We recommend commercially available nickel based solders. **Kera®NH** parts should not be soldered with gold or palladium solder. **Kera®NH** is also ideally suitable for laser-welding (commercial nickel based laser wire).

Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0,3 mm. It's recommended to sandblast the frames with minimum 110 µm of aluminium oxide with 3-4 bar and clean with steam cleaner.

Oxide firing is not mandatory but can be done as an option for 5 minutes at 980°C with vacuum (cleaning firing). The frame needs to be sandblasted with aluminium oxide 110 µm and 3-4 bar to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with NiCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

Single-use

The chemical and mechanical properties can only be guaranteed for new material.

Disposal Instructions

Consult the material safety data sheets or national regulations for disposal. Dispose of **Kera®NH** residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

Storage conditions

Temperature, humidity or light has no effect on the product properties.

Quantity

Please consider the label on the package.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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