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# Instruction for use Robur 400®

## Intended use

**Robur 400®** is a spring hard, beryllium- and nickel-free, cobalt based alloy for the manufacture of partial dentures. **Only for professional user!**

## General guidelines for handling

This instruction for use includes important processing steps and recommendations for **Robur 400®**.

## Indication

For the production of model cast works for removable restorations.

## Design

The design has to be made on an investment model (e.g. **Tener®Vest**, Eisenbacher Dentalwaren ED GmbH) under consideration of the standard rules of designs for partial dentures. The modelation should be done with wax that fire without leaving residues.

## Sprue design

The sprues should have a diameter of 3.5 mm with a length of 20-30 mm. Massive areas should additionally have a melt reservoir (ø 5mm). The cast hopper should be 4mm above the highest point of the partial denture wax design.

## Investing

Mixing ratios, expansion control, preheating processes are different depending on the investment material. Please refer to the processing guide of the corresponding investment material (e.g. **Tener®Vest**, Eisenbacher Dentalwaren ED GmbH). Recommended pre-heating temperature: 900 – 950°C.

## Melting and casting

**Robur 400®** should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material. Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and a so-called shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1490 °C.**

## Devesting

Let the muffle cool down to room temperature (ca. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 250 µm of aluminium oxide with 3-4 bar. Clean the partial denture base with a steam cleaner.

## Elaboration

The framework can be elaborated with standard carbide cutters or aluminium oxide stones and rubbers, look for smooth transitions.

## Soldering / Laser welding

**Robur 400®** can be soldered with all suitable solder. **Robur 400®** parts should not be soldered with gold or palladium solder. **Robur 400®** is also ideally suitable for laser-welding.

## Handling conditions / Safety

**Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 - DIN EN 149.**

## Contraindications and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the contra-indications and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorde d country.

## Single-use

The chemical and mechanical properties can only be guaranteed for new material.

## Disposal Instructions

Consult the material safety data sheets or national regulations for disposal. Dispose of **Robur 400®** residues and dust in an environmentally friendly manner. Grinding dust must not enter groundwater, water bodies or sewers. Address waste exchanges for recycling.

## Storage conditions

Temperature, humidity or light has no effect on the product properties.

## Quantity

Please consider the label on the package.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

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