Kera®S-Powder

DE - Gebrauchsanweisung Kera®S-Powder	2
EN - Instruction for use Kera®S-Powder	3
FR - Mode d'emploi Kera®S-Powder.....	3
CZ - Návod k použití Kera®S-Powder	4
EL - Οδηγίες χρήσης Kera®S-Powder	5
ES - Instrucciones de uso de Kera®S-Powder.....	6
HR - Upute za uporabu proizvoda Kera®S-Powder	7
HU - Használati útmutató a Kera®S-Powder termékhez	8
IT - Istruzioni per l'uso di Kera®S-Powder.....	9
PL - Instrukcja stosowania stopu Kera®S-Powder	10
PT - Instruções de utilização do Kera®S-Powder	11
RO - Instrucțiune de utilizare Kera®S-Powder	12
SE - Bruksanvisning för Kera®S-Powder.....	13
DA - Brugsanvisning til Kera®S-Powder	14
ET - Kasutusjuhend Kera®S-Powder	15
JA - ご使用方法 Kera®S-Powder	16
LT – Naudojimo instrukcija Kera®S-Powder.....	17
LV - Lietošanas pamācība Kera®S-Powder	18
NL - Gebruiksaanwijzing Kera®S-Powder	19
SK – Návod na použitie prášku Kera®S-Powder	20
SL – Navodila za uporabo Kera®S-Powder	21
TR - Kera®S-Powder Kullanım Talimatı.....	22

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
GBA / IFU Kera® S-Powder
Stand / Status 12/2024 as

0123

EN - Instruction for use Kera® S-Powder

Intended use
Kera®S-Powder is a medical device for the additive manufacturing of dental frameworks using the laser melting process.

Only for professional user (Dental Technician, Dentist).
The intended patient group provides for persons with partially or non-dentate jaws.

Indication
Crowns and bridges for metal-ceramics, frameworks for removable dentures, superstructures and bars.

Contraindication
In case of known allergic reactions to any of the ingredients.

CAD
The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm²). Sharp edges and undercuts should be avoided.

Stress-relief-heat-treatment
After the laser melting process, the building boards have to pass heat treatment to minimize stresses. For this purpose, a suitable furnace with inert gas (argon) or vacuum function should be used. Stress relieving without inert gas atmosphere can optionally be performed. Please consider that a treatment without inert gas atmosphere can lead to an increased oxide formation.

Description	Temperature [°C]	Time [min]
Stress-relief-heat-treatment with inert gas	850	30
Cooling 1 with inert gas	Cooling to 600, then open door	
Cooling 2, w/o inert gas	300 – room temperature	

Remove Parts from the build plate
After heat treatment and cooling the plate the restorations can be removed by band saw, rotary instruments or piers.

Soldering / Laser welding
Kera®S-Powder can be soldered with all suitable solder. Kera®S-Powder parts should not be soldered with gold or palladium solder. Kera®S-Powder is also ideally suitable for laser-welding.

Preparation before ceramic veneering
The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C with vacuum (cleaning firing). The frame needs to be sandblasted with aluminium oxide 110µm and 3-4 bar to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

Handling conditions / Safety
Metal powder or dust may cause irritation by inhalation and in contact with skin. During the handling with the Kera®S-Powder as well as while grinding and sandblasting of the produced units it's recommended to consider an adequate extraction system, goggles, gloves, protective clothing and a respirator with fine particle filter (type FFP3 – DIN EN 149). After working with the powder or restorations clean hands thoroughly.

Residual Risks and side effects
If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist about residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion
Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Disposal Instructions
Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions
Keep dry in a sealed container.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Eisenbacher Dentalwaren ED GmbH