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Stand / Status 05/2025 as

Instruction for use Kera®starSplint Bio

Mode d'emploi Kera®starSplint Bio

Intended use

Kera®starSplint Bio is a medical device for the milling of removable, occlusal dental splints.

Processing only by professional users! (dental technician, dentist).

Patient target group: Persons for therapeutic and orthodontic treatment with occlusal dental splints.

Material properties

Kera®starSplint Bio is characterised by excellent toughness, hydrolytic stability, heat resistance (max. 60°C) and chemical resistance.

Indication

- Patient-specific therapeutic occlusal splints to relieve the temporomandibular joint CMD
- Customised splints for bruxism

Contraindication

In case of intolerance to one of the ingredients.

Design

The virtual modelling is carried out with suitable CAD software, taking into account the rules of dental technology. The material thickness should not be less than 0.5 mm.

Milling

Please refer to the operating instructions and parameters of the respective CAM and milling machine manufacturer. The material is a thermoplastic and can be processed under water cooling. When milling without water cooling, an extraction unit must be used. Only use milling tools that are suitable for processing plastics and thermoplastics.

Cutting out the frameworks from the Blank

Cut off milled objects with suitable milling tools or cut-off wheels and grind supports.

Handling conditions / Safety

Dust may cause irritation by inhalation and in contact with skin. During grinding and sandblasting it's recommended to consider an adequate extraction system and goggles.

Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with **Kera®starSplint Bio** are extremely rare. In case of a proven allergy against an ingredient it should not be used for safety reasons. Please inform your dentist regarding the residual risk and side effects. Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country.

Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use

Used discs should be not processed for further use as a medical device.

Disposal Instructions

Please dispose residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions

The product must be stored in a dry place and protected from light.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Eisenbacher Dentalwaren ED GmbH