Robur 400®

DE - Gebrauchsanweisung Robur 400®	2
EN - Instruction for use Robur 400®	3
FR - Mode d'emploi Robur 400®	3
CZ - Návod k použití Robur 400®	4
EL - Οδηγίες χρήσης Robur 400®	5
ES - Instrucciones de uso de Robur 400®	6
HR - Upute za uporabu proizvoda Robur 400®	7
HU - Használati útmutató a Robur 400® termékhez	8
IT - Istruzioni per l'uso di Robur 400®	9
PL - Instrukcja stosowania stopu Robur 400®	10
PT - Instruções de utilização do Robur 400®	11
RO - Instrucțiune de utilizare Robur 400®	12
SE - Bruksanvisning för Robur 400®	13
DA - Brugsanvisning til Robur 400®	14
ET - Kasutusjuhend Robur 400®	15
JA - ご使用方法 Robur 400®	16
LT – Naudojimo instrukcija Robur 400®	17
LV - Lietošanas pamācība Robur 400®	18
NL - Gebruiksaanwijzing Robur 400®	19
SK – Návod na použitie zliatiny Robur 400®	20
SL – Navodila za uporabo Robur 400®	21
TR - Kullanım Talimatı Robur 400®	22

DE - Gebrauchsanweisung Robur 400®

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
EN - Instruction for use Robur 400®

Intended use

Robur 400® is a medical device for the manufacturing of cast partial dentures.

Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

Indication

For the production of model cast works for removable restorations.

Contraindication

In case of known allergic reactions to any of the ingredients.

Design

The design has to be made on an investment model under consideration of the standard rules of designs for partial dentures. The modelation should be done with wax that fire without leaving residues.

Sprue design

The sprues should have a diameter of 3.5 mm with a length of 20-30 mm. Massive areas should additional have a melt reservoir (ø 5mm). The cast hopper should be 4mm above the highest point of the partial denture wax design.

Melting and casting

Robur 400® should be melted in a ceramic crucible. Please do not use graphite crucibles and no flux! Avoid the overheating of the melt. Prevent multiple casts of melt bottoms. The chemical and mechanical properties can only be guaranteed for new material.

Melting with open flame (acetylene / oxygen) and inductive melting: Once the cylinders are melted and a so-called shadow falls across the molten metal, before the oxide skin begins to split, start the casting. **Maximum temperature for casting: 1490 °C.** The preheating temperature of the muffle is approx. 950°C.

Devesting

Let the muffle cool down to room temperature (ca. 20°C), do not quench with water. Put the cooled muffle into water to avoid dust generation during the devesting. Sandblast the surface with 250 µm of aluminium oxide with 3-4 bar. Clean the partial denture base with a steam cleaner.

Elaboration

The framework can be elaborated with standard carbide cutters or aluminium oxide stones and rubbers, look for smooth transitions.

Soldering / Laser welding

Robur 400® can be soldered with all suitable solder. **Robur 400®** parts should not be soldered with gold or palladium solder. **Robur 400®** is also ideally suitable for laser-welding.

Handling conditions / Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country.

Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use

The chemical and mechanical properties can only be guaranteed for new material.

Disposal Instructions

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions

Temperature, humidity or light has no effect on the product properties.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Eisenbacher Dentalwaren ED GmbH