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GBA/IFU Kera®Vest

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Instruction for use KERA[®]-VEST

Phosphate bonded investment powder for crowns and bridges

KERA[®]-VEST is a fine corn investment powder for universal applications, for casting alloys as well as pressable ceramic.

KERA[®]-VEST is approved for the speed heating and for the conventional heating process.

Storage- and processing temperature ca. 20-23 °C

Mixing ratio
100 g Powder: 24 - 26 ml liquid
160 g Powder: 39 - 41 ml liquid

Concentrate of mixing liquid

	Inlays, Onlays	Crowns and bridges (max. 3 units)	Bridges > 3 Units Secondary crowns
NPM alloy (CoCr, NiCr)	80 %*	80 - 85 %*	90 - 100 %*
PM (Gold alloy)	50 %*	55 - 60 %*	60 - 65 %*
Pressable ceramic	60 %*	60 %*	-

*The liquid concentration is only for orientation. The final liquid concentration have to find out by test casts.

Mixing Rinse the mixing bowl and dry with a clean towel. Add the required quantity of liquid into the mixing bowl. Add the appropriate quantity of KERA[®]-VEST powder. Mix well with a spatula (30 s.), then mixing under vacuum for 1 min, hold the vacuum for further 10 s after mixing.

Higher room temperatures shorten the processing time of the investment material!
Tip: Increase the mixing liquid by 1 ml.

Investing It's recommended doing the investing with silicone ring ((Use muffle fleece or roughen surfaces!). Once an investing with metal ring would be made it's mandatory to use sufficient layer (1-2) of fleece inlays. The fill in to the casting ring should be made with the lowest power of the vibrator. Once the ring has been filled up the vibrator have to switch off. Processing time ca. 5 min.

Speed heating Ca. 25 minutes after investing the casting form can be put into the hot furnace with a temperature of 850-900°C. Do not open the furnace door within the first 20 min.! Risk of injury!

Ring size	Time [min.]
1	35
3	45
6	60

Conventional heating With the conventional heating, the full setting expansion will be reached, you should therefore, reduce the liquid concentration by 5% - 10%.

Put the set up casting form into the cold oven. Please keep in mind to set a holding time at 290°C and 580°C for each 30-45 min (holding time depends on ring size). Hold the final temperature (850°C-900°C) also for 30 - 45 min.

General heat rate ca. 3 - 5 °C/min with linear oven control.

Safety instructions

Warning

- Do not open the furnace door during heat up.
- Investment powder contents silicon. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Use a Respirator, type FFP3-EN149. Avoid dust during opening and during fill in to the mixing bowl.
- Remove the dust at your working place only with wet towels. Avoid dust during deinvesting, place the cooled down muffle in water until they are thoroughly moistened. During sandblasting use an extraction.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.