# Kera®Ti5-Disc

DE - Gebrauchsanweisung Kera®Ti5-Disc .....	2
EN - Instruction for use Kera®Ti5-Disc.....	3
FR - Mode d'emploi Kera®Ti5-Disc .....	3
CZ - Návod k použití Kera®Ti5-Disc .....	4
EL - Οδηγίες χρήσης Kera®Ti5-Disc.....	5
ES - Instrucciones de uso de Kera®Ti5-Disc .....	6
HR - Upute za uporabu proizvoda Kera®Ti5-Disc.....	7
HU - Használati útmutató a Kera®Ti5-Disc termékhez .....	8
IT - Istruzioni per l'uso di Kera®Ti5-Disc .....	9
PL - Instrukcja stosowania stopu Kera®Ti5-Disc .....	10
PT - Instruções de utilização do Kera®Ti5-Disc .....	11
RO - Instrucțiune de utilizare Kera®Ti5-Disc.....	12
SE - Bruksanvisning för Kera®Ti5-Disc.....	13
DA - Brugsanvisning til Kera®Ti5-Disc .....	14
ET - Kasutusjuhend Kera®Ti5-Disc .....	15
JA - ご使用方法 Kera®Ti5-Disc .....	15
LT – Naudojimo instrukcija Kera®Ti5-Disc .....	16
LV - Lietošanas pamācība Kera®Ti5-Disc.....	17
NL - Gebruiksaanwijzing Kera®Ti5-Disc .....	18
SK – Návod na použitie zliatiny Kera®Ti5-Disc.....	19
SL – Navodila za uporabo Kera®Ti5-Disc.....	20
TR - Kera®Ti5-Disc Kullanım Talimatı .....	21

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
# EN - Instruction for use Kera®Ti5-Disc

## Intended use

Kera®Ti5-Disc is a medical device for the milling of crowns and bridges.

## Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

## Indication

Crowns and bridges for metal-ceramics up to 16 units (max. 4 units in the anterior region and max. 3 units in the posterior region), superstructures, bars.

## Contraindication

In case of known allergic reactions to any of the ingredients.

## Milling

Kera®Ti5-Disc is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and CNC milling machine.

## CAD

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6-9 mm<sup>2</sup>). Sharp edges and undercuts should be avoided.

## Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

## Laser welding

Sufficient and secure connections can only be done by laser welding. The use of Ti-welding wire is an option.

## Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with 110 µm of Aluminium oxide with 2-3 bar and clean with steam cleaner. Consider the firing process with max. 800°C, otherwise lattice transformation might occur.

## Handling conditions / Safety

**Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.**

## Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with Titan alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

## Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

## Single-use

Used discs should be not processed for further use as a medical device.

## Disposal Instructions

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

## Storage conditions

Temperature, humidity or light has no effect on the product properties.

*Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.*

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**ED GmbH, DIN EN ISO 13485 uyarınca
Eisenbacher Dentalwaren ED GmbH