Produktinformation / Product-Information

DE - Gebrauchsanweisung KERA® CAST 3D

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EN - Instruction for use KERA® CAST 3D

Intended use

KERA® CAST 3D is a light-curing resin for 3D printing of objects to be burnt out and pressed by means of the 3D printing technologies DLP or MSLA (LCD). KERA® CAST 3D can be used in 3D printers with 385 nm or 405 nm wavelength.

Contraindication

KERA® CAST 3D is not a medical device. It is only suitable for 3D printing of objects to be burnt out and pressed. Eisenbacher Dentalwaren is not liable for any damage that may result from incorrect use of the resin.

Processing instructions:

Preparation

- We recommend wearing personal protective equipment (e.g. gloves, safety goggles) when processing.
- Make sure that the 3D printer is in proper condition and that there is no contamination on the build platform, the material tray or the exposure unit. Observe the respective instructions for use for the 3D printer used.
- Shake the closed resin bottle properly (at least 2 minutes) before use.
- Then carefully pour KERA® CAST 3D into the material tray of the 3D printer.
- If there is still KERA® CAST 3D in the material tray, you should mix it thoroughly with a suitable plastic or silicone spatula for at least 2 minutes. Make sure that no material runs over the edge of the tray.
- Heating of the material by the 3D printer is not recommended.
- 3D printing
- The optional layer thicknesses of KERA® CAST 3D are 0.025 mm to 0.15 mm (recommended layer thickness: 0.05 mm).
- Select the product-specific parameter settings for KERA® CAST 3D for your printing process.
- After the end of the 3D printing process, let the printed objects hang on the build platform in your 3D printer for 5 minutes to drip off.
- Then scrape the printed object from the build platform using a suitable tool.
- Make sure that the printed objects are not exposed to direct sunlight before post-processing.

Cleaning

- Clean the printed objects with isopropanol (min. 90 %) in a suitable container. (Follow the safety instructions and disposal recommendations of the isopropanol manufacturer).
- Do not expose the printed objects to isopropanol for longer than 4 to 5 minutes.

Drying

- Dry the printed objects carefully with low pressure air or allow the objects to air dry for at least 7 minutes.

Post-exposure (mandatory)

2 x 400 flashes in Otofash G171

2 x 6 minutes in commercially available cure and wash equipment

Danger and safety information

Hazard statements: H302 Harmful if swallowed, H315 Causes skin irritation, H317 May cause allergic skin reactions, H318 Causes serious eye damage, H373 May cause damage to organs through prolonged or repeated exposure, H412 Harmful to aquatic organisms, with long lasting effects.

Precautionary statements: P101 If medical advice is needed, have product container or label at hand, P102 Keep out of reach of children, P280 Wear suitable protective clothing, gloves and eye/face protection, P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if possible. Continue to rinse, P310 Immediately call a POISON CENTRE/doctor, P501 Dispose of in accordance with local regulations.

Disposal advice

Dispose of liquid material bound in isopropanol or partially cured material in accordance with local / regional / national and international regulations.

Storage conditions

The material should be stored in a tightly closed container, in a place protected from light, at 18 °C to 28 °C. The liquid as well as the container should not be exposed to direct sunlight. The liquid as well as the container should not be exposed to sunlight and only briefly to daylight, otherwise premature curing may occur.

Batch number / expiry date

You will find this information on the packaging. Do not use the product after the expiry date.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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