Kera®Soft-Disc

DE - Gebrauchsanweisung Kera®Soft-Disc.....	2
EN - Instruction for use Kera®Soft-Disc	3
FR - Mode d'emploi Kera®Soft-Disc.....	3
CZ - Návod k použití Kera®Soft-Disc	4
EL - Οδηγίες χρήσης Kera®Soft-Disc	5
ES - Instrucciones de uso de Kera®Soft-Disc.....	6
HR - Upute za uporabu proizvoda Kera®Soft-Disc.....	7
HU - Használati útmutató a Kera®Soft-Disc termékhez.....	8
IT - Istruzioni per l'uso di Kera®Soft-Disc.....	9
PL - Instrukcja stosowania stopu Kera®Soft-Disc.....	10
PT - Instruções de utilização do Kera®Soft-Disc	11
RO - Instrucțiune de utilizare Kera®Soft-Disc	12
SE - Bruksanvisning för Kera®Soft-Disc.....	13
DA - Brugsanvisning til Kera®Soft-Disc.....	14
ET - Kasutusjuhend Kera®Soft-Disc	15
JA - ご使用方法 Kera®Soft-Disc	16
LT – Naudojimo instrukcija Kera®Soft-Disc.....	17
LV - Lietošanas pamācība Kera®Soft-Disc	18
NL - Gebruiksaanwijzing Kera®Soft-Disc	19
SK – Návod na použitie zliatiny Kera®Soft-Disc.....	20
SL – Navodila za uporabo Kera®Soft-Disc	21
TR - Kera®Soft-Disc Kullanım Talimatı.....	22

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
Stand / Status 12-2025 as

EN - Instruction for use Kera®Soft-Disc

Intended use
Kera®Soft-Disc is a medical device for the milling of crowns and bridges. The milled framework must be sintered under argon gas. **Only for professional user (Dental Technician, Dentist).**

Indication
Crowns and bridges for metal-ceramics

Contraindication and Patient group
In case of known allergic reactions to any of the ingredients.
The product is intended exclusively for use in adult patients. The use in vulnerable groups of people - in particular children, pregnant women and breastfeeding women - is contraindicated, as an increased sensitivity to potentially released metal ions cannot be ruled out in these groups.

Preparation recommendation	Connector cross section, anterior	Connector cross section, posterior	Wall thickness	Margin thickness
step or pronounced chamfer	6 mm²	9 mm²	0,4 – 0,5 mm	0,2 mm

- For restorations > 3 units, a sinter support analogous to zirconium oxide must be considered.
- The shrinkage factor is indicated on the blank and must be set in the CAM software parameters.
- Zirconium oxide or carbide cutters with the milling strategy for zirconium oxide or sinter metal are suitable for milling.
- The milling process is carried out dry under aspiration or with water cooling.
- With water cooling, dry the milled parts for 30 min at 120°C before sintering.

Information for sintering process

- Before sintering, the objects must be cleaned of milling dust with a brush.
- Make sure that the objects are completely embedded in the sinterball bed (approx. 1 - 2 mm distance from the margin).
- Bridge pontics, bar parts, etc. must not be covered with sinter beads.
- The objects must be placed in such a way that sinter shrinkage is not impeded.
- If several restorations are sintered together, they must not touch each other.
- After cooling to < 50°C, the finished parts can be removed from the furnace.

Programmable sintering furnace with argon		Amann Girrbach	
Heating rate	10°C / min	Heating rate	Standard
Final temperature	1280 °C	Final temperature	Standard
Holding time at final temperature	1 h	Holding time at final temperature	Standard
Cooling from end-temperature	Unregulated up to 300°C in argon atmosphere, from 300°C in air up to 50°C	Cooling from end-temperature	Standard

Soldering
Kera®Soft-Disc can be soldered with all suitable solder. **Kera®Soft-Disc** parts should not be soldered with gold or palladium solder.

Preparation before ceramic veneering
The frameworks are ground with cross-cut carbide cutters; ensure smooth transitions; avoid material overlaps. The minimum thickness of the caps should not be less than 0.3 mm. It is recommended that the frameworks are blasted with at least 110 µm aluminium oxide at 3-4 bar and cleaned (steam cleaner). An oxide firing can optionally be carried out at 980°C under vacuum with a holding time of one minute. **Kera®Soft-Disc** can be veneered with standard dental ceramics for CoCr. Please observe the process steps of the respective ceramic manufacturers.

Handling conditions / Safety
Metal dust is harmful to health. Always wear your personal protective equipment when handling the product. For milling, cleaning, finishing and sandblasting use suction (HEPA H filter) and respirator with FFP3 filter - DIN EN 149, gloves and goggles. Do not use compressed air for cleaning.

Residual risks and side effects
Incompatibilities with CoCr alloys are extremely rare if these instructions for use are observed. In case of a proven allergy to a component of this alloy, it must not be used for safety reasons. In exceptional cases, electrochemically induced local irritations are described. When using different alloy groups in the patient's mouth, galvanic effects can occur. Please inform your dentist about contraindications and side effects.
Any serious incident that involve the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion
Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Disposal Instructions
Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions
Temperature, humidity or light has no effect on the product properties

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

 **Eisenbacher Dentalwaren ED GmbH**
Dr.-Konrad-Wiegand-Str. 9 – 63939 Wörth am Main – GERMANY
Phone: +49 9372 94040 Fax: +49 9372 940429
E-Mail: info@eisenbacher.de Web: www.eisenbacher.de

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