Kera®Ti5-Powder

DE - Gebrauchsanweisung Kera®Ti5-Powder	2
EN - Instruction for use Kera®Ti5-Powder	3
FR - Mode d'emploi Kera®Ti5-Powder.....	3
CZ - Návod k použití Kera®Ti5-Powder	4
EL - Οδηγίες χρήσης Kera®Ti5-Powder	5
ES - Instrucciones de uso de Kera®Ti5-Powder	6
HR - Upute za uporabu proizvoda Kera®Ti5-Powder.....	7
HU - Használati útmutató a Kera®Ti5-Powder termékhez	8
IT - Istruzioni per l'uso di Kera®Ti5-Powder.....	9
PL - Instrukcja stosowania stopu Kera®Ti5-Powder	10
PT - Instruções de utilização do Kera®Ti5-Powder	11
RO - Instrucțiune de utilizare Kera®Ti5-Powder	12
SE - Bruksanvisning för Kera®Ti5-Powder.....	13
DA - Brugsanvisning til Kera®Ti5-Powder	14
ET - Kasutusjuhend Kera®Ti5-Powder	15
JA - ご使用方法 Kera®Ti5-Powder	16
LT – Naudojimo instrukcija Kera®Ti5-Powder.....	17
LV - Lietošanas pamācība Kera®Ti5-Powder	18
NL - Gebruiksaanwijzing Kera®Ti5-Powder.....	19
SK – Návod na použitie prášku Kera®Ti5-Powder	20
SL – Navodila za uporabo Kera®Ti5-Powder	21
TR - Kera®Ti5-Powder Kullanım Talimatı	22

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
Eisenbacher Dentalwaren ED GmbH

Dr.-Konrad-Wiegand-Str. 9 – 63939 Woerth am Main – GERMANY

Phone: +49 / 93 72 / 94 04 – 0 Fax: +49 / 93 72 / 94 04 – 29

E-Mail: info@eisenbacher.de

Web: www.eisenbacher.de

Abroad Link
Castellana Business Center
C/Paseo de la Castellana 40, 8ª Planta
Madrid 28046, Spain

GBA / IFU Kera® Ti5-Powder
Stand / Status 12/2024 as

EN - Instruction for use Kera®Ti5-Powder

Intended use

Kera®Ti5-Powder is a medical device for the additive manufacturing of dental frameworks using the laser melting process.

Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

Indication

Crowns and bridges for metal-ceramics, frameworks for removable dentures, superstructures and bars.

Contraindication

In case of known allergic reactions to any of the ingredients.

Design

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0,3 mm. Choose a sufficient connector dimension (6 - 9 mm²). Sharp edges and undercuts should be avoided.

Laser melting system

Before filling the laser melting system, check whether the system is designed for titanium powder processing. Observe the plant manufacturer's instructions for starting up the plant. After the laser melting process, the building plate must be cleaned from powder in compliance with the safety instructions.

Thermal treatment

The build plate with objects, which has been cleaned from powder, must be heat treated in order to minimise stresses in the material and to ensure the typical material properties. An electrically heated furnace with inert gas (argon) or a high vacuum furnace (min. 10⁻³ mbar) should be used for this purpose. Stress relief annealing without inert gas is not recommended due to the strong oxide formation. Due to the strong affinity of titanium to oxygen and hydrogen, the structure can be damaged.

Description	Temperature	Time
Stress-relief-heat-treatment with inert gas	700 °C	60 min
Cooling 1 with inert gas	Cooling to 500°C, then open door	
Cooling 2, w/o inert gas	300°C – room temperature	

Remove Parts from the build plate

After heat treatment and cooling the plate restorations can be removed by band saw or rotary instruments.

Laser welding

Laser welding with the same type of joining material can be used to produce mechanically high-strength and corrosion-resistant joints.

Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0,3 mm. It's recommended to sandblast the frames with 110 µm of Aluminium oxide with 2-3 bar and clean with steam cleaner. Consider the firing process with max. 800°C, otherwise lattice transformation might occur. Only ceramics for titanium materials should be selected for ceramic veneering.

Handling conditions / Safety instructions

Avoid dust formation! Handling the product (opening the packaging and transfilling) can lead to dust formation. Metal dusts are flammable. Keep away from all sources of fire. Due to the oxygen affinity of titanium, there is a danger of burning and exploding, especially on high temperatures. Keep metal fire extinguishing powder ready. Metal powder or dust can cause irritation when inhaled and in contact with the skin. When grinding and blasting the units, as well as when handling the powder, make sure there is sufficient extraction and wear personal protective equipment (PPE) (safety goggles, protective gloves, protective clothing as well as a respirator with fine particle filter (type FFP3 - DIN EN 149)). In case of eye contact, rinse the affected eye with water or saline solution for at least 15 min. Clean hands thoroughly after working with the powder or restorations. Observe the safety data sheet!

Handling conditions / Safety

Metal powder or dust may cause irritation by inhalation and in contact with skin. During the handling with the Kera®S-Powder as well as while grinding and sandblasting of the produced units it's recommended to consider an adequate extraction system, goggles, gloves, protective clothing and a respirator with fine particle filter (type FFP3 – DIN EN 149). After working with the powder or restorations clean hands thoroughly.

Residual Risks and side effects

If the instructions are observed during the production processes, incompatibilities with non-precious dental alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist about residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Disposal Instructions

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions

Keep dry in a sealed original container. Environmental temperatures have no influence on the product quality.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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