Kera[®]-Line

DE - Gebrauchsanweisung Kera [®] -Line	2
EN - Instruction for use Kera [®] -Line	3
FR - Mode d'emploi Kera [®] -Line	3
CZ - Návod k použití Kera [®] -Line	4
EL - Οδηγίες χρήσης Kera [®] -Line	5
ES - Instrucciones de uso de Kera [®] -Line	6
HR - Upute za uporabu proizvoda Kera [®] -Line	7
HU - Használati útmutató a Kera [®] -Line termékhez	8
IT - Istruzioni per l'uso di Kera [®] -Line	9
PL - Instrukcja stosowania stopu Kera [®] -Line	10
PT - Instruções de utilização do Kera [®] -Line	11
RO - Instrucțiune de utilizare Kera [®] -Line	12
SE - Bruksanvisning för Kera [®] -Line	13
DA – Brugsanvisning for Kera [®] -Line	14
ET – Kasutusjuhend Kera [®] -Line	15
JP - ご使用方法 Kera [®] -Line	16
LT – Naudojimo instrukcija Kera [®] -Line	17
LV – Lietošanas pamācība Kera [®] -Line	18
NL – Gebruiksaanwijzing Kera [®] -Line	19
SK – Návod na použitie zliatiny Kera [®] -Line	20
SL – Navodila za uporabo Kera [®] -Line	21
TR -Kera [®] -Line Kullanım Talimatı	22

DE - Gebrauchsanweisung Kera®-Line

PRODUKTNAME PRODUCT NAME / NOM DU PRODUIT	Kera®-Line
BEZEICHNUNG DESCRIPTION / DESCRIPTIF	NEM Dental-Fräslegierung auf Kobaltbasis für Metall-Keramik, Typ 4 / NPM Cobalt-based dental milling alloy for metal-ceramic, type 4 / Alliage de fraisage dentaire NEM à base de cobalt pour métal-céramique, type 4
ABMESSUNG DIMENSION / DIMENSIONS	Ø 12,0 mm x 1000 mm

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0123

EN - Instruction for use Kera®-Line

Intended use
Kera®-Line is a medical device for the milling of crowns and bridges.

Only for professional user (Dental Technician, Dentist).

Indication
Crowns and bridges for metal-ceramics, superstructures.

Contraindication and Patient group
In case of known allergic reactions to any of the ingredients.

The product is intended exclusively for use in adult patients. The use in vulnerable groups of people - in particular children, pregnant women and breastfeeding women - is contraindicated, as an increased sensitivity to potentially released metal ions cannot be ruled out in these groups.

Milling
Kera®-Line is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine.

CAD
The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm²). Sharp edges and undercuts should be avoided.

Cutting out the frameworks from the milling alloy
Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Soldering / Laser welding
Kera®-Line can be soldered with all suitable solder. **Kera®-Line** parts should not be soldered with gold or palladium solder. **Kera®-Line** is also ideally suitable for laser-welding.

Preparation before ceramic veneering
The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980°C (cleaning firing). The frame needs to be sandblasted again to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder, please consider the instruction for use of the manufacturer.

Handling conditions / Safety
Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Residual risks and side effects
If the instructions are observed during the production processes, incompatibilities with CoCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion
Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use
Used milling alloys should be not processed for further use as a medical device.

Disposal Instructions
Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions
Temperature, humidity or light has no effect on the product properties

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Eisenbacher Dentalwaren ED GmbH