de	en	fr	es	it	nl	da	no	sv	fi	pt	cs	pl
-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------	-----------------

KERA®-TRAY

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Quantitative Angaben:
In 100 g there are 19 g acrylates/methacrylates and 80 g Glas enthalten.

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Times:
Working at daylight 20 minutes
Polymerisation in an ultraviolet or LED light-curing unit, depending on the unit in use. The top of the tray, on the model 3-5 minutes
The underside of the tray, off the model 3-5 minutes
Trimming:
- Trim the crown with tungsten carbide burrs and smooth with sandpaper.

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KERA®-TRAY

Instructions for use

Patient target group
All persons for whom individual trays are required for the treatment of individual cases.

Intended user
Dentist

Composition Qualitative data:
Oligomeric urethane acrylates, vinyl ester resin, 2,2-dimethoxy-1,2-diphenylethan-1-one or by 2,2-dimethoxy-1,2-diphenylethan-1-one with 1,2,3-Dioxol-1,7-trimethyl-bicyclo-(2,2)-heptan-2-yl-2,6-dimethyl-2-phenylacetophenon, (+)-2,3-Dioxol-1,7-trimethyl-bicyclo-(2,2)-heptan-2-yl-2,6-dimethyl-2-phenylacetophenon, 2-hydroxy-2-methyl-1-phenylpropan-1-one, 2-ethylhexyl-4-(dimethylamino)benzoate, 2-ethylhexyl-4-(dimethylamino)benzoate, 2,2-diethoxyacetophenone, amorphous silica, dyes (Pigment Red 146 = CI 12485).

Quantitative data:
In 100 g there are 19 g acrylates/methacrylates and 80 g glass.

Intended purpose
KERA®-TRAY is a light-curing tray material in the form of prefabricated brackets to make rigid, custom impression trays and bite-blocks, as well as for other indications. KERA®-TRAY can be cured in all ultraviolet and halogen light-curing units. The baseplates are supplied in a re-usable light-proof box that can be used for transporting or storing models, etc.

Indication
- Custom impression tray
- Bite-block
- Custom for setting up dentures
- Base for supporting pin retention
- Support for prosthetic or telescoping crowns
- Control base for denture groups

Contraindications
The material contains methacrylates. Do not use in patients with an allergy to methacrylates or other substances contained in the product.

Performance characteristics of the product
KERA®-TRAY polymerises with light to form a solid and compact polymer.

Special precautions
Avoid any contact between the eyes and skin with uncured material. Use the system with the inhibited layer.

Always use a dust extractor when trimming the tray or when using the material.

Preparation/Manufacture
- Mark the outline of the future tray on the model.
- Block out the undercuts on the model.
- To create a spacer for the impression material, use a wax or a spacer.
- If the patient is very dry, soak the model briefly in water or apply eye isolation or an alginate isolation.

- Remove the baseplate from the box and grind the model.
- Once the model has been blocked out, attach the baseplate to it and cut it to shape using a dust extractor.
- Use the remaining material to form a handle using a wax or the tray.
- Place the model and tray in a light-curing unit and cure without vacuum.
- After the tray has been cured, remove the vacuum polymerisation tray from the underside. If the tray has been cured without a vacuum, remove the tray from the underside.

Times:
Working at daylight 20 minutes
Polymerisation in an ultraviolet or LED light-curing unit, depending on the unit in use. The top of the tray, on the model 3-5 minutes
The underside of the tray, off the model 3-5 minutes

Trimming:
- Trim the crown with tungsten carbide burrs and smooth with sandpaper.

Hygiene:
Immerse impressions taken with KERA®-TRAY customs trays, in a disinfectant solution based on glutaraldehyde. The time depends on the disinfectant, manufacturer's directions.

The finished medical product is used by the dentist, e.g. for individual trays with impression material.

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