Gebrauchsanweisung TENER VEST

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GBA/IFU TENER VEST

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Instruction for Use TENER VEST

Phosphate-bonded investment material for model casting.

TENER VEST is approved for the "speed heating" and for the conventional heating process.

Storage- and processing temperature	ca. 20-23 °C
Mixing ratio	Silicone-duplication: 100 g powder : 20-22 ml fluid → 400 g : 80 ml
	Gel-duplication: 100 g powder : 18-20 ml fluid → 400 g : 72 ml

Concentrate of mixing liquid

	Partial denture with brackets	Combination prostheses (depends on size and extension)	Secondary elements, one piece casting
Liquid concentration	65 - 75 %*	75 - 85 %*	85 - 100 %*

*The liquid concentration is only for orientation. The final liquid concentration must be determined by means of test casts.

Mixing

Rinse the mixing bowl and dry with a clean towel. Add the required quantity of liquid into the mixing bowl. Add the appropriate quantity of **TENER VEST** powder. Mix well with a spatula (for 20 s), then mix under vacuum for 1 min (about 350 rpm), hold the vacuum for further 10 s after mixing.

Higher room temperatures shorten the processing time of the investment material!
Tip: Increase the mixing liquid by 1 ml.

Investing

The fill-in into the casting ring should be made at the lowest power vibrator. Once the casting ring has been filled up, the vibrator must be switched off. Processing time approx. 5 min, setting time approx. 9 min.

The overbedding can be mixed with distilled water.

Speed heating

Ca. 25 min after fill-in, the casting ring can be put into the pre-heated furnace at a temperature of 850 - 900°C. Hold final temperature of 900-950°C for at least 45 min. Do not open the furnace door within the first 20 minutes! Risk of injury!

Conventional heating

With the conventional heating, the full setting expansion will be reached, you should therefore reduce the liquid concentration by 5% - 10%.

Put the set-up casting form into the cold oven. Please keep in mind to set a holding time at 280°C and 580°C of 45 - 60 min (depending on the casting ring size and number of rings). Hold the final temperature of 900°C – 950°C for at least 45 min.

General heating rate approx. 3 - 7 °C/min with linear oven control.

Safety notes

Warning

- Do not open the furnace door during heat-up.
- Investment powder contents silicone. Do not inhale the dust! Risk of lung damage (silicosis, lung cancer). Recommendation: Use a Respirator, type FFP3-EN149. Avoid dust during opening and during fill-in to the mixing bowl.
- Remove the dust at your working place only with wet towels.
- Avoid dust during deinvesting, place the cooled-down muffle in water until it is thoroughly moistened.
- Use an extraction system during sandblasting.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.